




EQUINE ARTHROSCOPY



R RITA
LEIBINGER
MEDICAL

Curved Hook Knife

Design by Dr. Hans Wilderjans

163-200-20



www.LEIBINGER-MEDICAL.com

Indication:

Chronic tenosynovitis of the digital flexor tendon sheath with thickening and constriction of the palmar annular ligament.

Goal:

Section of a thickened and constricted annular ligament using a tenoscopic approach


This knife is specially designed to facilitate easy cutting of the palmar annular ligament in cases of chronic tenosynovitis of the digital flexor tendon sheath in the horse.



Specifications

- Curved shaft and bigger hook to facilitate good grip and deep cut in thickened ligaments
- New handle will facilitate accurate placement in hand
- Sharp blade ensures easy cutting in firm chronic thickened ligament

Arthroscopy Rongeurs



Schlessinger Rongeur

3x10mm serrated jaw, 5"shaft, straight


163-404-755



Schlessinger Rongeur

4x10mm serrated jaw, 5"shaft, straight


163-404-756



Wilde Ethmoid Rongeur

4mm jaw, 5,5"shaft, straight

163-404-768



McIlwraith Peapod Rongeur

2.5x6mm, 5,5"shaft, straight

163-404-782



McIlwraith Peapod Rongeur

up-angle, Medium Jaw, 5,5" Shaft

163-404-784

Arthroscopy Rongeurs



Ferris-Smith Rongeur

straight, 3x10mm jaw, 7"

163-404-791



Ferris-Smith Rongeur

straight, 4x10mm jaw, 7"

163-404-792



Ferris-Smith Rongeur

straight, 6x10mm bite, 7"

163-404-794



Ferris-Smith Rongeur

straight, 6x12mm bite, 7"

163-404-795



Ferris-Smith Rongeur

up-angle, 4x10mm bite, 7"

163-404-802



Ferris-Smith Rongeur

down-angle, 4x10mm bite, 7"

163-404-812

Arthroscopy Currettes



Brun's Curette 9"

oval, straight, size 00

163-320-831

oval, curved, size 00

163-320-851



Brun's Curette 9"

oval, straight, size 0

163-320-832

oval, curved, size 0

163-320-852



Brun's Curette 9"

oval, straight, size 1

163-320-833

oval, curved, size 1

163-320-853



Brun's Curette 9"

oval, straight, size 2

163-320-834

oval, curved, size 2

163-320-854



Brun's Curette 9"

oval, str., size 000

163-320-870



Brun's Curette 9"

oval, str., size 00

163-320-871



Brun's Curette 9"

oval, str., size 0

163-320-872

Arthroscopy Knives and Hooks



Knife

45°, 23 cm

163-1271-126



Scalpel Knife

23 cm

163-1271-128



Button-end Knife

23 cm

163-1271-132



Dandy Nerv Hook

1 mm, 19 cm

163-1408-00

Hook Probe, Bone Pick and Elevators



Mcllwraith Hook Probe

heavy, 27 cm

163-1271-130



Equine Bone Pick

25 cm

163-1271-325



Mcllwraith Periosteal Elevator

19 cm, 5mm tip

163-1404-870



Foerner Knife Elevator

16 cm, 5mm tip, cvd.

35/3714-875