

LAPAROSCOPY



Laparoscopy

Designed by RITA LEIBINGER



1. Length

Length and Diameter were adjusted for veterinary medicine. This enables ergonomic working.



2. Rotatable



Rotatable working ends enable ergonomic and accurate working.

3. Exchangeable working ends

The handle is compatible to different working ends. This enables economical working.



4. Cleaning

The working ends can be cleaned very easily. All instruments are detachable. For sterilizing purposes the instruments may be detached into 2-parts (\varnothing 3 mm instruments) or 3-parts (\varnothing 3,5 / 5 / 10 mm instruments) of components. All items are compatible with each other and can be sterilized easily.



PEEK Handpiece

with ratchet or free mechanism
with HF connection, rotatable

167-4001-00



Shaft only

with screw plug
Ø5mm, 250mm
for 167-4502-25, 167-4503-25 and 167-4504-25

167-4505-25

with screw plug
Ø5mm, 310mm
for 167-4502-31, 167-4503-31 and 167-4504-31

167-4505-31



Maryland Grasping Forceps

serrated
Ø3mm, 250mm

167-4301-25

Ø3mm, 310mm

167-4301-31



Grasping Forceps

fenestrated, serrated
Ø3mm, 250mm

167-4302-25

Ø3mm, 310mm

167-4302-31



Biopsy Forceps Spoon

oval
Ø3mm, 250mm

167-4303-25

Ø3mm, 310mm

167-4303-31



Metzenbaum Scissors

insulated, curved
double action, Ø3mm, 250mm

167-4304-25

Ø3mm, 310mm

167-4304-31



Metzenbaum Scissors

insulated, curved
Ø5mm, 250mm, insert only
double action

167-4502-25

Ø5mm, 310mm

167-4502-31



Grasping Forceps

pointed, with spoon
serrated,
Ø5mm, 250mm, insert only

167-4503-25

Ø5mm, 310mm

167-4503-31



Biopsy Forceps Spoon

oval
Ø5mm, 250mm
insert only

167-4504-25

Ø5mm, 310mm

167-4504-31



Grasping Forceps

serrated
Ø5mm, 250mm
insert only

167-4610-25

Ø5mm, 310mm

167-4610-31



Screwshaft Trocar Set, PEEK, complete

trocar blunt and trocar sharp,
Length 60mm

167-0030-50



Endoscopes

Endoscopes

rigid, autoclavable 132°/273°F,
opening angle 0°, stainless steel



Product Code	Diameter	Length
161-5005-25	5 mm	250 mm
161-5005-30	5 mm	300 mm
161-5010-25	10 mm	250 mm
161-5010-30	10 mm	300 mm