


INSTRUMENTS FOR MICRO SURGERY



Rita Leibinger Medical is proud to present the latest development (advancement) in orthopedic instrument technology. With our continuous development of specialty lines of instruments we want to introduce to you this new line of Small Bone Orthopedic Instruments.

The **Small Bone Orthopedic Instruments** are offered in 72 various Instruments, and will offer surgeons the choice to create their own SBOI Set. We offer the SBOI Sterilization Container for 5 – 10 and 20 Instruments.

Original Size (15 cm)





[1] **Bone Lever**
cvd. up



[2] **Bone Lever**
cvd. side



[3] **Bone Hook**
sharp 1 prong



[4] **Bone Hook**
blunt 1 prong



[5] **Langenbeck**



[6] **Perio. Elevator**
slight cvd.



[7] **Bone Curette**



[8] **Bone Curette**



[9] **Gouge**



[10] **Bone File**
cvd., coarse



[11] **Bone File**
cvd., fine



[12] **Bone Awl**



[13] **Bone Pa-
cker**
ball tip



[14] **Joseph
Hook**
2 prong



[15] **Senn Miller
Hook**
sharp, 3 prong



[16] **Senn Miller
Hook**
blunt, 3 prong



[17] **Senn Miller
Hook**
sharp, 2 prong



[18] **Meyerding
Retractor**



[19] **Meyerding
Retractor**



[20] **Meyerding
Retractor**
concave



[21] **Scalpel
Handle**
angled



[22] **Scalpel
Handle**
straight



[23] **Cushing
Retractor**



[24] **Osteotome**



[25] **Hohmann
Retractor**

Instrument cases



[26] Instrument
Cassette #5
for 5 instruments



[27] Instrument
Cassette #10
for 10 instruments



[28] Instrument Cassette #20
for 20 instruments

N°	Code	Description
[1]	170-1038-03	Bone Lever cvd. up, 3 mm tip
[2]	170-1037-03	Bone Lever, cvd. side, 3 mm tip
[3]	170-1170-01	Bone Hook, sharp, 1 prong
[4]	170-1225-07	Bone Hook, blunt, 1 prong
[5]	170-1235-06	Langenbeck Retractor, 10x6 mm
	170-1236-04	Langenbeck Retractor, 11x4 mm
[6]	170-1025-02	Periosteal Elevator, slight cvd., 2 mm tip
	170-1026-03	Periosteal Elevator, slight cvd., 3 mm tip
	170-1027-04	Periosteal Elevator, slight cvd., 4 mm tip
	170-1028-05	Periosteal Elevator, slight cvd., 5 mm tip
	170-1029-06	Periosteal Elevator, slight cvd., 6 mm tip
	170-1030-07	Periosteal Elevator, slight cvd., 7 mm tip
	170-1031-08	Periosteal Elevator, slight cvd., 8 mm tip
[7] & [8]	170-1040-02	Bone Curette, Cup Size 2 mm
	170-1045-03	Bone Curette, Cup Size 3 mm
	170-1050-04	Bone Curette, Cup Size 4 mm
	170-1055-05	Bone Curette, Cup Size 5 mm
[9]	170-1082-02	Gouge, 2 mm tip
	170-1083-03	Gouge, 3 mm tip
	170-1084-04	Gouge, 4 mm tip
	170-1085-05	Gouge, 5 mm tip
	170-1086-06	Gouge, 6 mm tip
	170-1087-07	Gouge, 7 mm tip
	170-1088-08	Gouge, 8 mm tip
	170-1090-10	Gouge, 10 mm tip
	170-1092-12	Gouge, 12 mm tip
[10]	170-1104-04	Bone File, 4 mm, cvd., coarse
	170-1105-06	Bone File, 6 mm, cvd., coarse
[11]	170-1124-04	Bone File, 4 mm, cvd., fine
	170-1125-06	Bone File, 6 mm, cvd., fine
[12]	170-1135-05	Bone Awl, 5 mm tip
	170-1137-25	Bone Awl, 2.5 mm tip
[13]	170-1155-02	Bone Packer, 2 mm ball tip
	170-1156-03	Bone Packer, 3 mm ball tip
	170-1157-04	Bone Packer, 4 mm ball tip
	170-1158-05	Bone Packer, 5 mm ball tip
[14]	170-1210-03	Joseph Hook, sharp, 2-prong, 3 mm wide
	170-1211-02	Joseph Hook, sharp, 2-prong, 2 mm wide
	170-1212-07	Joseph Hook, sharp, 2-prong, 7 mm wide
	170-1215-06	Joseph Hook, sharp, 2-prong, 4 mm wide
	170-1216-04	Joseph Hook, bl., 2-prong, 4 mm wide
	170-1218-05	Joseph Hook, sh., 2-prong, 8 mm wide
[15]	170-1220-08	Senn Miller Hook, sh., 3-prong, 8mm cvd., 11mm wide
[16]	170-1222-08	Senn Miller Hook, bl., 3-prong, 8mm cvd., 11mm wide
[17]	170-1223-08	Senn Miller Hook, sh., 2-prong, 8mm cvd., 7mm wide
	170-1226-08	Senn Miller Hook, bl., 2-prong, 8mm cvd., 7mm wide
[18]	170-1232-04	Meyering Retractor, 6x4 mm
[19]	170-1233-04	Meyering Retractor, 6x4 mm, concave
[20]	170-1234-04	Meyering Retractor, 10x4 mm
[21]	170-1305-03	Scalpel Handle, No. 3, angled
[22]	170-1300-03	Scalpel Handle, No. 3, straight
[23]	170-1266-06	Cushing Retractor, 6 mm, rigid
	170-1268-08	Cushing Retractor, 8 mm, rigid, concave
	170-1269-08	Cushing Retractor, 8 mm, rigid, flat
	170-1270-10	Cushing Retractor, 10 mm, rigid
	170-1273-12	Cushing Retractor, 12 mm, rigid
[24]	170-1062-02	Osteotome, 2 mm tip
	170-1063-03	Osteotome, 3 mm tip
	170-1064-04	Osteotome, 4 mm tip
	170-1065-05	Osteotome, 5 mm tip
	170-1066-06	Osteotome, 6 mm tip
	170-1067-07	Osteotome, 7 mm tip
	170-1068-08	Osteotome, 8 mm tip
	170-1070-10	Osteotome, 10 mm tip
	170-1072-12	Osteotome, 12 mm tip
[25]	170-1005-06	Hohmann Retractor, 6 mm tip
	170-1010-09	Hohmann Retractor, 9 mm tip
[26]	150-0400-05	Instrument Cassette for five hand instruments
[27]	150-0400-10	Instrument Cassette for ten hand instruments
[28]	150-0400-20	Instrument Cassette for twenty hand instruments