


IBIAL TUBEROSITY ADVANCEMENT



TTA ESY Sets




TTA ESY Instruments BOX

without contents

BOX112

(similar picture)



TTA ESY Implants BOX

without contents

BOX116



**TTA Instruments
& Implants BOXES**

without contents

BOX115

Cages

Titanium



Product Code	Size (mm)
230-2162-313	3/13
230-2162-316	3/16
230-2162-319	3/19
230-2162-4512	4.5/12
230-2162-4515	4.5/15
230-2162-4518	4.5/18
230-2162-616	6/16
230-2162-619	6/19
230-2162-622	6/22
230-2162-7516	7.5/16
230-2162-7519	7.5/19
230-2162-7522	7.5/22
230-2162-919	9/19
230-2162-922	9/22
230-2162-925	9/25
230-2162-10519	10.5/19
230-2162-10522	10.5/22
230-2162-10525	10.5/25
230-2162-1222	12/22
230-2162-1225	12/25
230-2162-1228	12/28
230-2162-1525	15/25
230-2162-1528	15/28
230-2162-1531	15/31

ESY TTA Plates

Titanium



Product Code	Specification	
	Holes	Style
230-155-03	2	3
230-155-04	2	4
230-155-05	2	5
230-155-06	3	6
230-155-07	3	7
230-155-08	3	8

Forks

Titanium



Prongs	Product Code
3	230-2161-03
4	230-2161-04
5	230-2161-05
6	230-2161-06
7	230-2161-07
8	230-2161-08

Plates Standard



Holes	Product Code
3	230-2160-03
4	230-2160-04
5	230-2160-05
6	230-2160-06
7	230-2160-07
8	230-2160-08

Fork Inserter

246-002-00



Self Tapping Screws

2.4 mm Hex Head,
Titanium



Measurements	
Head dia	4.0 mm
Thread dia	2.4 mm
Core dia	1.7 mm
Thread pitch	1.0 mm

Product Code	Length (mm)
245-224-06	6
245-224-08	8
245-224-10	10
245-224-12	12
245-224-14	14
245-224-16	16
245-224-18	18
245-224-20	20
245-224-22	22
245-224-24	24
245-224-26	26
245-224-28	28
245-224-30	30
245-224-32	32
245-224-34	34
245-224-36	36
245-224-38	38
245-224-40	40

Self Tapping Screws

2.7 mm Hex Head
Titanium



Measurements	
Head dia	5.0 mm
Thread dia	2.7 mm
Core dia	1.9 mm
Thread pitch	1.0 mm

Product Code	Length (mm)
245-227-08	8
245-227-10	10
245-227-12	12
245-227-14	14
245-227-16	16
245-227-18	18
245-227-20	20
245-227-22	22
245-227-24	24
245-227-26	26
245-227-28	28
245-227-30	30
245-227-32	32
245-227-34	34
245-227-36	36
245-227-38	38
245-227-40	40

Self Tapping Screws

3.5 mm Hex Head
Titanium



Hex Head	
Head dia	6.0 mm
Thread dia	3.5 mm
Core dia	2.4 mm
Thread pitch	1.25 mm

Product Code	Length (mm)
245-235-10	10
245-235-12	12
245-235-14	14
245-235-16	16
245-235-18	18
245-235-20	20
245-235-22	22
245-235-24	24
245-235-26	26
245-235-28	28
245-235-30	30
245-235-32	32
245-235-34	34
245-235-36	36
245-235-38	38
245-235-40	40

TTA Spacers

Titanium



Product Code	Size (mm)
246-008-02	2
246-008-04	4
246-008-06	6

Plate Bender

246-003-00



T-Handle
246-004-00



Spreaders



Width (mm)	Product Code
------------	--------------

6	246-005-06
---	-------------------

9	246-005-09
---	-------------------

12	246-005-12
----	-------------------

Drill Guides for Forks



Holes	Product Code
-------	--------------

4	246-001-04
---	-------------------

8	246-001-08
---	-------------------

Pin



Size	Product Code
------	--------------

1.9	246-006-19
-----	-------------------

Drill Sleeves



Product Code	Size (mm)
164-0070-24	1.8/2.0
164-0070-27	2.7/2.0
164-0070-35	3.5/2.5

Drill Bits



Product Code	Length (mm)	Spec.
148-0080-18	1.8	straight shank
148-0080-20	2.0	straight shank
148-0080-25	2.5	straight shank
148-0081-18	1.8	AO QC
148-0081-20	2.0	AO QC
148-0081-25	2.5	AO QC

Depth Gauge

for 2.4 | 2.7 | 3.5 screws
maximal measurement 40 mm

164-2427-35



Screwdriver Handle

w. Quick Coupling
sterilizable up to 134 °C

128-0910-24



Screwdriver Hex Insert

128-0900-20	2.4 mm
128-0900-25	2.7/3.5 mm



Holding Sleeve for Screwdriver

128-0940-24	2.4 mm
128-0940-25	2.7/3.5 mm



Plate Holding Forceps

16 cm angled

164-0050-16



TTA RAPID® Spreader

with tensioning and fixation screw, 16 cm

132-4080-16



Boneholding Forceps

21,5 cm - with spin lock

128-0525-21

