



MUSTANG



MUSTANG Battery Driven Drill

Handpiece, for large bones, without battery

125-0140-00



MUSTANG Jacob Chuck

with key
max. Ø 6.0 mm drill

125-0141-00



MUSTANG Battery only

NiMh Battery not autoclavable

125-0142-00



MUSTANG Charger for one Battery

125-0143-00



Battery Steril Funnel

for TERRIER and MUSTANG handpiece
autoclavable

125-0117-00



non-sterile battery
non-sterile battery must be
inserted into the sterile
handpiece.

Therefore you have a ...

sterile battery funnel
through which you insert the
non-sterile battery into the

sterile handpiece