




# TIBIAL TUBEROSITY ADVANCEMENT



## TTA ESY Sets

### TTA ESY Implants BOX

*without contents*

**BOX116**



### Sterilization Container

310x190x130mm  
blue (image)

**150-5401-30**

green

**150-5402-30**



Cages

Titanium



Product Code	Size (mm)
230-2162-313	3/13
230-2162-316	3/16
230-2162-319	3/19
230-2162-4512	4.5/12
230-2162-4515	4.5/15
230-2162-4518	4.5/18
230-2162-616	6/16
230-2162-619	6/19
230-2162-622	6/22
230-2162-7516	7.5/16
230-2162-7519	7.5/19
230-2162-7522	7.5/22
230-2162-919	9/19
230-2162-922	9/22
230-2162-925	9/25
230-2162-10519	10.5/19
230-2162-10522	10.5/22
230-2162-10525	10.5/25
230-2162-1222	12/22
230-2162-1225	12/25
230-2162-1228	12/28
230-2162-1525	15/25
230-2162-1528	15/28
230-2162-1531	15/31

ESY TTA Plates

Titanium



Product Code	Specification	
	Holes	Style
230-155-03	2	3
230-155-04	2	4
230-155-05	2	5
230-155-06	3	6
230-155-07	3	7
230-155-08	3	8

Forks

Titanium



Prongs	Product Code
3	230-2161-03
4	230-2161-04
5	230-2161-05
6	230-2161-06
7	230-2161-07
8	230-2161-08

Plates Standard



Holes	Product Code
3	230-2160-03
4	230-2160-04
5	230-2160-05
6	230-2160-06
7	230-2160-07
8	230-2160-08

Fork Inserter

246-002-00



Plate Bender

246-003-00



### Self Tapping Screws

2.4 mm Hex Head,  
Titanium



Measurements	
Head dia	4.0 mm
Thread dia	2.4 mm
Core dia	1.7 mm
Thread pitch	1.0 mm

Product Code	Length (mm)
245-224-06	6
245-224-08	8
245-224-10	10
245-224-12	12
245-224-14	14
245-224-16	16
245-224-18	18
245-224-20	20
245-224-22	22
245-224-24	24
245-224-26	26
245-224-28	28
245-224-30	30
245-224-32	32
245-224-34	34
245-224-36	36
245-224-38	38
245-224-40	40

### Self Tapping Screws

2.7 mm Hex Head  
Titanium



Measurements	
Head dia	5.0 mm
Thread dia	2.7 mm
Core dia	1.9 mm
Thread pitch	1.0 mm

Product Code	Length (mm)
245-227-08	8
245-227-10	10
245-227-12	12
245-227-14	14
245-227-16	16
245-227-18	18
245-227-20	20
245-227-22	22
245-227-24	24
245-227-26	26
245-227-28	28
245-227-30	30
245-227-32	32
245-227-34	34
245-227-36	36
245-227-38	38
245-227-40	40

### Self Tapping Screws

3.5 mm Hex Head  
Titanium



Hex Head	
Head dia	6.0 mm
Thread dia	3.5 mm
Core dia	2.4 mm
Thread pitch	1.25 mm

Product Code	Length (mm)
245-235-10	10
245-235-12	12
245-235-14	14
245-235-16	16
245-235-18	18
245-235-20	20
245-235-22	22
245-235-24	24
245-235-26	26
245-235-28	28
245-235-30	30
245-235-32	32
245-235-34	34
245-235-36	36
245-235-38	38
245-235-40	40

### TTA Spacers

Titanium



Product Code	Size (mm)
246-008-02	2
246-008-04	4
246-008-06	6

### Drill Guides for Forks



Holes	Product Code
4	246-001-04
8	246-001-08

### Pin



Size	Product Code
1.9	246-006-19

### Drill Sleeves



Product Code	Size (mm)
164-0070-24	1.8/2.0
164-0070-27	2.7/2.0
164-0070-35	3.5/2.5

### Drill Bits



Product Code	Ø (mm)	Spec.
148-0080-18	1.8	straight shank
148-0080-20	2.0	straight shank
148-0080-25	2.5	straight shank
148-0081-18	1.8	AO QC
148-0081-20	2.0	AO QC
148-0081-25	2.5	AO QC

### Depth Gauge

for 2.4 | 2.7 | 3.5 screws  
maximal measurement 40 mm

**164-2427-35**



### Screwdriver Handle

w. Quick Coupling  
sterilizable up to 134 °C

**128-0910-24**



**Screwdriver Hex Insert**

<b>128-0900-20</b>	2.4 mm
<b>128-0900-25</b>	2.7/3.5 mm



**Holding Sleeve for Screwdriver**

<b>128-0940-24</b>	2.4 mm
<b>128-0940-25</b>	2.7/3.5 mm



**Plate Holding Forceps**

16 cm angled

**164-0050-16**



**TTA RAPID® Spreader**

with tensioning and fixation screw, 16 cm

**132-4080-16**



**Boneholding Forceps**

21,5 cm - with spin lock

**128-0525-21**

