



LeiLOX



TPLO Locking System



R **RITA**
LEIBINGER
MEDICAL

TPLO Set

TPLO Set
142-0100-00



Product Code	Description	Quantity
142-0100-10	TPLO Tray Implantats	1,00
142-0100-20	TPLO Tray Instruments	1,00
128-2735-10	Star-Drive 10 Screwdriver shaft	1,00
164-2700-00	Locking Drill Guide 2.7 Stainless Steel	2,00
164-3500-00	Locking Drill Guide 3.5 Stainless Steel	2,00
144-1010-10	K-Draht 1.0, 10 cm single Trokar, round end	2,00
164-2735-60	Depth Gauge	1,00
148-0080-20	Twist Drill Bone round shaft, 2.0 mm	1,00
148-0081-20	Twist Drill Bone AO shaft, 2.0 mm	1,00
148-0080-25	Twist Drill Bone round shaft, 2.5 mm	1,00
148-0081-25	Twist Drill Bone AO shaft, 2.5 mm	1,00
128-2720-27	Compressional Drill Guide, neutral/load, 2.7 mm, for 2,7 mm screws and 2.0 mm drills	1,00
128-2535-35	Compressional Drill Guide, neutral/load, 3.5 mm, for 3.5 mm screws and 2.5 mm drills	1,00
128-0900-25	Screwdriver Insert, for 2,7/3,5 Screw, Star-Drive	1,00
128-0940-25	Holding sleeve for 2,7/3,5mm Screws	1,00
247-0102-00	Screwdriver Handle Silicone, straight	1,00
142-1127-00	LeiLox TPLO Locking Plate, 2.7 mm, right, 46 mm, Stainless Steel	2,00
142-1127-10	LeiLox TPLO Locking Plate, 2.7 mm, left, 46 mm, Stainless Steel	2,00
142-1135-00	LeiLox TPLO Locking Plate, 3.5 mm, right, 66 mm, Stainless Steel	2,00
142-1135-10	LeiLox TPLO Locking Plate, 3.5 mm, left, 66 mm, Stainless Steel	2,00
142-1135-50	LeiLox TPLO Locking Plate, 3.5 mm breit, right, 75 mm, Stainless Steel	2,00
142-1135-60	LeiLox TPLO Locking Plate, 3.5 mm breit, left, 75 mm, Stainless Steel	2,00
245-627-16	Bone Screw Cortex, 2.7, 16 mm length, self-tapping, Star-Drive, stainless st.	3,00
245-627-18	Bone Screw Cortex, 2.7, 18 mm length, self-tapping, Star-Drive, stainless st.	3,00
245-627-20	Bone Screw Cortex, 2.7, 20 mm length, self-tapping, Star-Drive, stainless st.	3,00
245-627-22	Bone Screw Cortex, 2.7, 22 mm length, self-tapping, Star-Drive, stainless st.	3,00
245-627-24	Bone Screw Cortex, 2.7, 24 mm length, self-tapping, Star-Drive, stainless st.	3,00
245-627-26	Bone Screw Cortex, 2.7, 26 mm length, self-tapping, Star-Drive, stainless st.	3,00
245-627-28	Bone Screw Cortex, 2.7, 28 mm length, self-tapping, Star-Drive, stainless st.	3,00
245-627-30	Bone Screw Cortex, 2.7, 30 mm length, self-tapping, Star-Drive, stainless steel	3,00
245-627-32	Bone Screw Cortex, 2.7, 32 mm length, self-tapping, Star-Drive, stainless steel	3,00
245-627-34	Bone Screw Cortex, 2.7, 34 mm length, self-tapping, Star-Drive, stainless steel	3,00
245-635-16	Bone Screw Cortex, 3.5, 16 mm length, self-tapping, Star-Drive, stainless steel	3,00
245-635-18	Bone Screw Cortex, 3.5, 18 mm length, self-tapping, Star-Drive, stainless steel	3,00
245-635-20	Bone Screw Cortex, 3.5, 20 mm length, self-tapping, Star-Drive, stainless steel	3,00
245-635-22	Bone Screw Cortex, 3.5, 22 mm length, self-tapping, Star-Drive, stainless steel	3,00
245-635-24	Bone Screw Cortex, 3.5, 24 mm length, self-tapping, Star-Drive, stainless steel	3,00
245-635-26	Bone Screw Cortex, 3.5, 26 mm length, self-tapping, Star-Drive, stainless steel	3,00
245-635-28	Bone Screw Cortex, 3.5, 28 mm length, self-tapping, Star-Drive, stainless steel	3,00
245-635-30	Bone Screw Cortex, 3.5, 30 mm length, self-tapping, Star-Drive, stainless steel	3,00
245-635-32	Bone Screw Cortex, 3.5, 32 mm length, self-tapping, Star-Drive, stainless steel	3,00
245-635-34	Bone Screw Cortex, 3.5, 34 mm length, self-tapping, Star-Drive, stainless steel	3,00
242-127-16	Locking Bone Screw 2.7 mm, 16 mm lenght, Star-Drive, self-tapping, stainless steel	5,00
242-127-18	Locking Bone Screw 2.7 mm, 18 mm lenght, Star-Drive, self-tapping, stainless steel	5,00
242-127-20	Locking Bone Screw 2.7 mm, 20 mm lenght, Star-Drive, self-tapping, stainless steel	5,00

242-127-22	Locking Bone Screw 2.7 mm, 22 mm length, Star-Drive, self-tapping, stainless steel	5,00
242-127-24	Locking Bone Screw 2.7 mm, 24 mm length, Star-Drive, self-tapping, stainless steel	5,00
242-127-26	Locking Bone Screw 2.7 mm, 26 mm length, Star-Drive, self-tapping, stainless steel	5,00
242-127-28	Locking Bone Screw 2.7 mm, 28 mm length, Star-Drive, self-tapping, stainless steel	5,00
242-127-30	Locking Bone Screw 2.7 mm, 30 mm length, Star-Drive, self-tapping, stainless steel	5,00
242-127-32	Locking Bone Screw 2.7 mm, 32 mm length, Star-Drive, self-tapping, stainless steel	5,00
242-127-34	Locking Bone Screw 2.7 mm, 34 mm length, Star-Drive, self-tapping, stainless steel	5,00
242-127-36	Locking Bone Screw 2.7 mm, 36 mm length, Star-Drive, self-tapping, stainless steel	5,00
242-127-38	Locking Bone Screw 2.7 mm, 38 mm length, Star-Drive, self-tapping, stainless steel	5,00
242-127-40	Locking Bone Screw 2.7 mm, 40 mm length, Star-Drive, self-tapping, stainless steel	5,00
242-127-42	Locking Bone Screw 2.7 mm, 42 mm length, Star-Drive, self-tapping, stainless steel	5,00
242-127-44	Locking Bone Screw 2.7 mm, 44 mm length, Star-Drive, self-tapping, stainless steel	5,00
242-127-46	Locking Bone Screw 2.7 mm, 46 mm length, Star-Drive, self-tapping, stainless steel	5,00
242-135-16	Locking Bone Screw 3.5 mm, 16 mm length, Star-Drive, self-tapping, stainless steel	5,00
242-135-18	Locking Bone Screw 3.5 mm, 18 mm length, Star-Drive, self-tapping, stainless steel	5,00
242-135-20	Locking Bone Screw 3.5 mm, 20 mm length, Star-Drive, self-tapping, stainless steel	5,00
242-135-22	Locking Bone Screw 3.5 mm, 22 mm length, Star-Drive, self-tapping, stainless steel	5,00
242-135-24	Locking Bone Screw 3.5 mm, 24 mm length, Star-Drive, self-tapping, stainless steel	5,00
242-135-26	Locking Bone Screw 3.5 mm, 26 mm length, Star-Drive, self-tapping, stainless steel	5,00
242-135-28	Locking Bone Screw 3.5 mm, 28 mm length, Star-Drive, self-tapping, stainless steel	5,00
242-135-30	Locking Bone Screw 3.5 mm, 30 mm length, Star-Drive, self-tapping, stainless steel	5,00
242-135-32	Locking Bone Screw 3.5 mm, 32 mm length, Star-Drive, self-tapping, stainless steel	5,00
242-135-34	Locking Bone Screw 3.5 mm, 34 mm length, Star-Drive, self-tapping, stainless steel	5,00
242-135-36	Locking Bone Screw 3.5 mm, 36 mm length, Star-Drive, self-tapping, stainless steel	5,00
242-135-38	Locking Bone Screw 3.5 mm, 38 mm length, Star-Drive, self-tapping, stainless steel	5,00
242-135-40	Locking Bone Screw 3.5 mm, 40 mm length, Star-Drive, self-tapping, stainless steel	5,00
242-135-42	Locking Bone Screw 3.5 mm, 42 mm length, Star-Drive, self-tapping, stainless steel	5,00
242-135-44	Locking Bone Screw 3.5 mm, 44 mm length, Star-Drive, self-tapping, stainless steel	5,00
242-135-46	Locking Bone Screw 3.5 mm, 46 mm length, Star-Drive, self-tapping, stainless steel	5,00



TPLO Locking Plates

LeiLox TPLO Locking Plate



LeiLox TPLO Locking Plate, 3.5 mm,
left, 75 mm, stainless steel

142-1135-60



LeiLox TPLO Locking Plate, 3.5 mm,
right, 75 mm, stainless steel

142-1135-50



LeiLox TPLO Locking Plate, 3.5 mm,
left, 66 mm, stainless steel

142-1135-10



LeiLox TPLO Locking Plate, 3.5 mm,
rechts, 66 mm, stainless steel

142-1135-00



LeiLox TPLO Locking Plate, 2.7 mm,
left, 46 mm, stainless steel

142-1127-10



LeiLox TPLO Locking Plate, 2.7 mm,
right, 46 mm, stainless steel

142-1127-00

Screws (Star-Drive,

Locking Bone Screw Star-Drive

Stainless Steel
for LeiLox Locking Plates
self-tapping, self-holding



Length (mm)	Ø (mm)	Prod. Code	Ø (mm)	Prod. Code
16	2.7	242-127-16	3.5	242-135-16
18	2.7	242-127-18	3.5	242-135-18
20	2.7	242-127-20	3.5	242-135-20
22	2.7	242-127-22	3.5	242-135-22
24	2.7	242-127-24	3.5	242-135-24
26	2.7	242-127-26	3.5	242-135-26
28	2.7	242-127-28	3.5	242-135-28
30	2.7	242-127-30	3.5	242-135-30
32	2.7	242-127-32	3.5	242-135-32
34	2.7	242-127-34	3.5	242-135-34
36	2.7	242-127-36	3.5	242-135-36
38	2.7	242-127-38	3.5	242-135-38
40	2.7	242-127-40	3.5	242-135-40
42	2.7	242-127-42	3.5	242-135-42
44	2.7	242-127-44	3.5	242-135-44
46	2.7	242-127-46	3.5	242-135-46

Bone Screw Star-Drive (Not Locking)

stainless steel,
self-tapping, self-holding



Length (mm)	Ø (mm)	Prod. Code	Ø (mm)	Prod. Code
16	2.7	245-627-16	3.5	245-635-16
18	2.7	245-627-18	3.5	245-635-18
20	2.7	245-627-20	3.5	245-635-20
22	2.7	245-627-22	3.5	245-635-22
24	2.7	245-627-24	3.5	245-635-24
26	2.7	245-627-26	3.5	245-635-26
28	2.7	245-627-28	3.5	245-635-28
30	2.7	245-627-30	3.5	245-635-30
32	2.7	245-627-32	3.5	245-635-32
34	2.7	245-627-34	3.5	245-635-34

Screwdriver Shaft Star-Drive 10

self-holding

128-2735-10



Screwdriver Handle

Silicone

sterilizable up to 134°C

247-0102-00



Torque Limiting Screwdriver

Torque 3 Nm

dark blue, Silicone, straight, AO
sterilizable up to 134°C

247-0104-00



K-Wire

Single Trokar



1.0 mm x 100 mm **144-1010-10**

Bone Drills



Ø (mm)	Connection	Product Code
2.0	round shaft	148-0080-20
2.5	round shaft	148-0080-25
2.0	AO shaft	148-0081-20
2.5	AO shaft	148-0081-25

Depth Gauge

164-2735-60



Compressional Drill Guide

neutral/load

128-2720-27

2.7 mm, for screws Ø 2.7 mm and drills Ø 2.0 mm



128-2535-35

3.5 mm, for screws Ø 3.5 mm and drills Ø 2,5mm



Locking Drill Guide

stainless steel



for 2.7 mm screws
164-2700-00



for 3.5 mm screws
164-3500-00